

Prepared Media Plates for Microbiological Testing

Easy Plate™ EC (Medi•Ca EC)

Inclusivity and Exclusivity

● *E. Coli*

According to the research by DNP

Strain name	Source	Easy Plate EC (Medi•Ca EC)	Chromogenic agar media	BGLB broth
<i>Escherichia coli</i>	NBRC102203	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC12062	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC12433	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC12734	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC13500	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC15034	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC3972	+	(Blue-Purple)	+
<i>Escherichia coli</i>	ATCC25922	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC3301	+	(Blue-Purple)	+
<i>Escherichia coli</i>	NBRC13540	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC3543	+	(Blue-Purple)	+
<i>Escherichia coli</i>	NBRC3806	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC3991	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC13898	+	(Navy-Blue)	+
<i>Escherichia coli</i>	NBRC14129	+	(Blue-Purple)	+
<i>Escherichia coli</i>	D0099	+	(Navy-Blue)	+
<i>Escherichia coli</i>	D0100	+	(Navy-Blue)	+
<i>Escherichia coli</i>	D0101	+	(Navy-Blue)	+
<i>Escherichia coli</i>	D0102	+	(Navy-Blue)	+
<i>Escherichia coli</i>	D0104	+	(Navy-Blue)	+
<i>Escherichia coli</i> O157:H7※	ATCC43895	-	(Red-Purple)	
<i>Escherichia coli</i> O26:NM	RIMD05091876	+	(Blue-Purple)	

※ (+): Positive, (-): Negative, () Color of colony

※ *Escherichia coli* O157 is detected as non-*E. coli* coliform and forms red-purple colonies because it does not produce β-glucuronidase.

● Non-*E. coli* Coliform

According to the research by DNP

Strain name	Source	Easy Plate EC (Medi•Ca EC)	Chromogenic agar media	BGLB broth
<i>Citrobacter amalonaticus</i>	NBRC13547	+	(Red-Purple)	+
<i>Citrobacter freundii</i>	NBRC12681	+	(Red-Purple)	+
<i>Citrobacter koseri</i>	NBRC105690	+	(Red-Purple)	+
<i>Cronobacter sakazakii</i>	D0003	+	(Red-Purple)	+
<i>Buttiauxella noackiae</i>	D0077	+	(Pink)	+
<i>Enterobacter asburiae</i>	D0029	+	(Red-Purple)	+
<i>Enterobacter aerogenes</i>	NBRC13534	+	(Red-Purple)	+
<i>Enterobacter amnigenus</i>	D0037	+	(Red-Purple)	+
<i>Enterobacter cloacae</i>	D0030	+	(Red-Purple)	+
<i>Enterobacter cloacae</i>	ATCC222	+	(Red-Purple)	+
<i>Enterobacter cloacae</i>	D0033	+	(Pink)	+
<i>Escherichia blattae</i> (<i>Shimwellia blattae</i>)※	NBRC105725	-	(Pink)	-
<i>Escherichia fergusonii</i>	NBRC102419	+	(Red-Purple)	+
<i>Escherichia hermannii</i>	NBRC105704	+	(Red-Purple)	+
<i>Escherichia vulneris</i>	NBRC102420	+	(Red-Purple)	+
<i>Klebsiella oxytoca</i>	D0032	+	(Pink)	+
<i>Klebsiella oxytoca</i>	NBRC105695	+	(Red-Purple)	+
<i>Klebsiella pneumoniae</i>	ATCC13883	+	(Red-Purple)	+
<i>Kluyvera cryocrescens</i>	NBRC102467	+	(Red-Purple)	+
<i>Leclercia adecarboxylata</i>	NBRC102595	+	(Red-Purple)	+
<i>Pantoea agglomerans</i>	D0004	+	(Pink)	-
<i>Rahnella aquatilis</i>	D0053	+	(Red-Purple)	-
<i>Rahnella sp.</i>	D0023	+	(Red-Purple)	-
<i>Raoultella terrigena</i>	D0022	+	(Red-Purple)	+

※ (+): Positive, (-): Negative, () Color of colony

※ *Escherichia blattae* cannot be detected because it does not produce β-galactosidase.

※ Numbers starting with D indicate strains that were isolated from foods by Dai Nippon Printing Co., Ltd.

● Non-Coliform

According to the research by DNP

Strain name	Source	Easy Plate EC (Medi・Ca EC)	Chromogenic agar media	BGLB broth
<i>Achromobacter xylosoxidans</i>	NBRC15126	-	- (White)	-
<i>Bacillus amylolycheferiens</i>	D0015	-	-	-
<i>Bacillus cereus</i>	NBRC3836	-	-	-
<i>Bacillus coagulans</i>	NBRC12583	-	-	-
<i>Bacillus licheniformis</i>	D0010	-	-	-
<i>Bacillus megaterium</i>	NBRC15308	-	-	-
<i>Bacillus subtilis</i>	D0021	-	-	-
<i>Bacillus thuringiensis</i>	NBRC3951	-	-	-
<i>Corynebacterium variabile</i>	NBRC15286	-	-	-
<i>Edwardsiella tarda</i>	NBRC105688	-	- (White)	-
<i>Enterococcus faecalis</i>	ATCC29212	-	-	-
<i>Enterococcus faecium</i>	NBRC100486	-	-	-
<i>Kocuria rhizophila</i>	D0008	-	-	-
<i>Lactobacillus casei</i>	D0025	-	-	-
<i>Lactobacillus delbrueckii</i>	NBRC3202	-	-	-
<i>Lactococcus lactis</i>	D0026	-	-	-
<i>Leuconostoc mesenteroides</i>	D0057	-	-	-
<i>Micrococcus luteus</i>	NBRC3333	-	-	-
<i>Proteus hauseri</i>	NBRC3851	-	- (White)	-
<i>Proteus mirabilis</i>	NBRC105697	-	- (White)	-
<i>Pseudomonas aeruginosa</i>	NBRC3899	-	- (White)	-
<i>Pseudomonas mendocina</i>	NBRC14162	-	- (White)	-
<i>Pseudomonas sp.</i>	D0054	-	- (White)	-
<i>Salmonella enterica</i>	NBRC105726	-	- (White)	+
<i>Serratia liquefaciens</i>	D0027	-	- (White)	-
<i>Serratia marcescens</i> ※	NBRC102204	+ (Red-Purple)	+ (Pink)	-
<i>Serratia rubidaea</i> ※	NBRC12973	+ (Red-Purple)	+ (Pink)	-
<i>Staphylococcus epidermidis</i>	NBRC100911	-	-	-
<i>Staphylococcus aureus</i>	D0072	-	-	-
<i>Staphylococcus aureus</i>	ATCC25923	-	-	-
<i>Staphylococcus carnosus</i>	D0086	-	-	-
<i>Staphylococcus gallinarum</i>	D0061	-	-	-
<i>Staphylococcus intermedius</i>	ATCC29663	-	-	-
<i>Staphylococcus saprophyticus</i>	D0009	-	-	-
<i>Staphylococcus simulans</i>	NBRC109714	-	-	-
<i>Staphylococcus sp.</i>	D0058	-	-	-
<i>Staphylococcus xylosus</i>	NBRC109770	-	-	-
<i>Yersinia frederiksenii</i>	D0052	-	- (White)	-

※ (+): Positive, (-): Negative, () Color of colony

※ Some *Serratia* species may be classified as coliform because of the β -galactosidase production.

※ Numbers starting with D indicate strains that were isolated from foods by Dai Nippon Printing Co., Ltd.

※ These results are just examples and may vary depending on the condition of the strain, culture conditions, etc.

※ Incubation Condition

Easy Plate EC (Medi・Ca EC) : 35±1°C, 24 hours

Chromogenic agar media : 35±1°C, 24 hours

BGLB broth : 35±1°C, 24-48 hours

Dilution solution: Butterfield's phosphate-buffered diluent)



Kikkoman Biochemifa Company

2-1-1 Nishi-Shinbashi, Minato-ku, Tokyo 105-0003, Japan

E-mail: biochemifa@mail.kikkoman.co.jp

URL: <https://biochemifa.kikkoman.com/e/>

※ The information contained in this catalog may change without notice.

©2021 Kikkoman Corp. (6012W211201)