

CleanTrust™

Surface Protein Allergen Test Swab



The CleanTrust™ Surface Protein Allergen Test Swab is a ready to use rapid test swab with high sensitivity, used to detect protein residues on surface, and thereby determine if cleaning has been effective in eliminating proteins, including potentially allergenic proteins. After smearing the surface of the object, if protein is present, the solution will change from green to purple. The faster and darker the color change, the higher the protein content. It is important to note that a residue may contain a mixture of allergenic and non-allergenic proteins, as well as non-protein materials.


Benefits

- Visual interpretation, no instrument needed
- High sensitivity, detects as little as 5 µg of protein
- Suitable for protein-based allergens monitoring

Sensitivity



- 5 µg after 15 minutes incubation at 37°C

Procedure



Result Interpretation




 Visual Interpretation



Order Information



Product No.	Product Name	Specification
YL0200	CleanTrust™ Surface Protein Allergen Test Swab	100 pcs / box

 Telephone: (02) 9540 2288 or 1300 628 104
 E-mail: info@fmcgis.com.au
 Address: Unit 11A, 1 – 3 Endeavour Road
 Caringbah NSW 2229 AUSTRALIA