

The CleanTrust Surface Protein Test Swab is a ready to use rapid test swab used to detect protein residues on cleaned surfaces. After smearing the surface of the object, if protein is present, the solution will change from green to purple. The faster and darker the color change, the higher the protein content.


Benefits

- Easy to use, no extra reagent needed
- · Visual interpretation, no instrument needed
- Quick and stable, result in 1-10 minutes

Sensitivity


- Room temperature (20-25°C):
- 80µg after 1 minute
- 50µg after 5 minutes
- 20µg after 10 minutes


Procedure





Result Interpretation





Order Information



Product No.	Product Name	Specification
YL0100	CleanTrust [™] Surface Protein Test Swab	100 pcs / box

- © Telephone: (02) 9540 2288 or 1300 628 104
- Address: Unit 11A, 1 3 Endeavour Road Caringbah NSW 2229 AUSTRALIA